

### iNtouch Latex Micro Textured Sterile Surgical Gloves

INTW10

## ouch®

### **PRODUCT INFORMATION**

ltem Number	INTW10
Product Colour	Natural
Size Range	5.5-9
Fingertip Thickness	0.23 mm
Palm Thickness	0.20 mm
Typical Weight	12.70 g
Shelf Life	5 Years

### FEATURES

Sterilised: Sterilised by irradiation method

Hand Specific: Hand Specific reduces hand fatigue

Micro Roughened: Provides excellent feel and instrument grip in both wet and dry conditions

**Comfort:** Manufactured on ergonomically designed formers from an advanced synthetic formulation, offers ultimate comfort and tactile sensitivity

Cuff: Long Beaded Cuff

Manufactured From: Natural rubber latex, which can cause allergic reactions, including anaphylactic responses.

### SECTORS

Medical, Vet and Dental

### **STANDARDS**

### AQL 0.65

This product is classified as Class IIa Medical Device in compliance with the requirements of Medical Device Regulation (EU) 2017/745 & UK MDR 2002 (SI2002 No. 618 as amended) (UKCA)

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425



Compliant with European Standards EN455: Medical Gloves

MD

CLASS II: MEDICAL DEVICE

# Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-5: 2016











### CAUTION

**STANDARDS** 

for Single Use parts 1,2,3&4

EN ISO 374-4 & EN ISO 374-5

STERILE R

Contains latex which may cause allergic reactions, including anaphylactic response.

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individual

Polyco Healthline, South Fen Business Park, South Fen Road, Bourne, Lincolnshire, PE10 0DN Tel: 0333 320 8550 • E: sales@polycohealthline.com • www.polycohealthline.com

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555



PACKING

50 gloves per dispenser

4 dispensers per case

48 cases per pallet

INTW10

## iNtouch

### STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. You are advised to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear.

### **DISPOSAL INSTRUCTIONS**

Dispose of responsibly and in line with local guidance.

### **INNER PACKAGING**



Size	Inner Packaging Barcode
5.5	9555256608049
6	9555256608056
6.5	9555256608063
7	9555256608070
7.5	9555256608087
8	9555256608094
8.5	9555256608100
9	9555256608117

Size	Outer Packaging Barcode
5.5	9555256607967
6	9555256607974
6.5	9555256607981
7	9555256607998
7.5	9555256608001
8	09555256608018
8.5	9555256608025
9	9555256608032

#### PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thicknesses is approximate and based on an agreed manufacturing tolerance. All thickness and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555