

GN31 Handsafe Natural Latex Powder Free Examination



PRODUCT INFORMATION

ltem Number	GN31
Product Colour	Natural
Size Range	XS-XL
Fingertip Thickness	0.11 mm
Palm Thickness	0.10 mm
Typical Weight	5.50 g
Shelf Life	3 Years

FEATURES

Chemical Resistant: Achieves type C chemical resistance according to EN374-1, for chemical protection against at least 1 chemical

Virus/Non Virus Standard: Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

Chlorinated: Chlorinated to reduce latex protein content

Flexible material: Flexible material for improved feel and fit reducing fatigue during extended wear

Manufactured from: Manufactured from natural rubber latex

Powder Free: Reduced risk of dust contamination **Textured Finish:** Textured finish for improved grip

SECTORS

Hair and Beauty • Medical, Vet and Dental

STANDARDS

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5

AQL 1.5

This product is classified as Class I Medical Device in compliance with the requirements of Medical Device Regulation (EU) 2017/745 & UK MDR 2002 (SI2002 No.618 as amended) (UKCA)



STANDARDS

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4













2777

VIRUS





CAUTION

Contains latex which may cause allergic reactions, including anaphylactic response.

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555



GN31 Handsafe Natural Latex Powder Free Examination Glove



STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. You are advised to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear.

DISPOSAL INSTRUCTIONS

Dispose of responsibly and in line with local guidance.

PACKING

100 gloves per dispenser 10 dispensers per case 48 cases per pallet

INNER PACKAGING





OUTER PACKAGING





Size	Inner Packaging Barcode
Small	5024951023548
Medium	5024951023555
Large	5024951023562
XL	5024951023579

Inner Packaging		
Inner dimensions (LxWxH)	230 x 120 x 65mm	

Size	Outer Packaging Barcode
Small	5024951023593
Medium	5024951023609
Large	5024951023616
XL	5024951023623

Outer Packaging		
Outer Dimensions (LxWxH)	335 x 250 x 238mm	

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thicknesses is approximate and based on an agreed manufacturing tolerance. All thickness and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

 $For further technical information, please contact Customer Services on 0333\,320\,8550\, or the \, Technical \, team \, on \, 0333\,320\,8555\, or \, the \, Technical \, team \, on \, 0333\,320\, or \, the \, Technical \, team \, on \, 0333\,320\, or \, the \, Technical \, team \, on \, 0333\,320\, or \, the \, Technical \, team \, on \, 0333\,320\, or \, the \, Technical \, team \, on \, 0333\,320\, or \, the \, Technical \, team \, on \, 0333\, or \, the \, Technical \, team \, on \, 0333\, or \, the \, Technical \, team \, on \, 0333\, or \, the \, Technical \, team \, on \, 0333\, or \, the \, Technical \, team \, on \, 0333\, or \, the \, Technical \, team \, on \, 0333\, or \, the \, Technical \, team \, on \, 0333\, or \, the \, Technical \, team \, on \, 0333\, or \, the \, Technical \, team \, on \, 0333\, or \, the \, Technical \, team \, on \, the \, Technical \, team \, on \, t$