

Finite P Indigo AF Nitrile Accelerator Free Powder Free Exam Glove



PRODUCT INFORMATION

ltem Number	MFNP100
Product Colour	Purple
Size Range	XS - XL
Fingertip Thickness	0.08 mm
Palm Thickness	0.06 mm
Typical Weight	4.20 g
Shelf Life	3 Years

FEATURES

Chemical Resistant: Achieves type C chemical resistance according to EN ISO 374-1, for chemical protection against at least 1 chemical

Virus/Non Virus Standard: Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

Fit: Unique shape delivers a glove with improved fit

Accelerator Free: Accelerator free removing the risk of contact dermatitis caused by these chemicals

Comfort: Nitrile technology provides an exceptionally soft and pliable material

Latex Free: Latex Free reducing the risk of hypersensitivity and allergic reactions

Powder Free: Reduced risk of dust contamination

Sensitive Skin: No latex proteins - eliminates protein sensitisation

Suitability: Suitable for use with cytotoxic/chemotherapy

SECTORS

Automotive • Construction & Utility • Engineering and Manufacturing • Food Manufacturing and Preparation • Janitorial • Oil, Gas & Renewables



STANDARDS

AQL 1.5

This product is classified as Class I Medical Device in compliance with the requirements of Medical Device Regulation (EU) 2017/745 & UK MDR 2002 (SI2002 No.618 as amended) (UKCA)

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Tested in compliance with Commission Regulation (EU) No 10/2011 for Food contact

Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4

Tested to European standards EN420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5











VIRUS



2777



For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555



Finite P Indigo AF Nitrile Accelerator Free Powder Free **Exam Glove**



STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. You are advised to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear.

DISPOSAL INSTRUCTIONS

Dispose of responsibly and in line with local guidance.

PACKING

100 gloves per dispenser 10 dispensers per case 60 cases per pallet

INNER PACKAGING





Inner Packaging Barcode

<u>21</u>
PAP

XS	5010699508973
Small	5010699508980
Medium	5010699508997
Large	5010699509000
XL	5010699509017

Inner Packaging		
Inner dimensions (LxWxH)	230 x 120 x 68mm	

OUTER PACKAGING





Size	Outer Packaging Barcode
XS	05010699816115
Small	05010699816122
Medium	05010699816139
Large	05010699816146
XL	05010699816153

Outer Packaging		
Outer Dimensions (LxWxH)	352 x 240 x 248mm	

PLEASE NOTE:

Size

Weight does not form part of the specification and is subject to change. Palm thicknesses is approximate and based on an agreed manufacturing tolerance. All thickness and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555